


Bulk

Single-use suction tubes for Otology

Our sterile bulk packs contain 2 bags of 25 suction tubes that can be ordered in various lengths and diameters. This packaging saves valuable time and is more environmentally friendly.

- Luer lock connection
- Angled cannula made of stainless steel
- Large selection of lengths and diameters

| Art.-No. | Length (mm) | Diameter (mm) | Color |
|----------|-------------|---------------|-------|
| 302010U | 70 | 1.0 | cream |
| 302012U | 70 | 1.2 | pink |
| 302014U | 70 | 1.4 | red |
| 302016U | 70 | 1.6 | white |
| 302020U | 70 | 2.0 | green |
| 303014U | 80 | 1.4 | red |
| 303020U | 80 | 2.0 | green |
| 303025U | 80 | 2.5 | blue |
| 303030U | 80 | 3.0 | blue |
| 303040U | 80 | 4.0 | blue |
| 304014U | 100 | 1.4 | red |
| 304020U | 100 | 2.0 | green |
| 304025U | 100 | 2.5 | blue |
| 304030U | 100 | 3.0 | blue |



2 x 25 pieces per box, sterile packed

Not suitable for use in a sterile environment!



