

Black Line

Ultimate Edition



Quality you can rely on




Quality you can rely on

Black Line

Ultimate Edition

Benefit from the advantages of a Diamond-Like Carbon coating.
Excellent wear resistance and highest precision in your OR.



- Low wear thanks to DLC coating
- Low-reflective surface
- Good biocompatibility
- Many Wavecut models available

Wavecut is a special blade design which ensures an optimal dissection of soft tissue



Diamond-Like Carbon, DLC

The term DLC stands for a group of amorphous, carbon-based diamond-like layers.

DLC coatings are gaining in importance thanks to their outstanding wear, friction and corrosion-reducing properties.

DLC coatings are used where extreme hardness, low wear and minimal friction are required.



Black Line

Ultimate Edition

Scissors

- 10-532-12B Tungsten carbide, straight, pointed, 11.5 cm
- 10-532-12WB Tungsten carbide, straight, pointed, 11.5 cm, Wavecut
- 10-535-12B Tungsten carbide, curved, pointed, 11.5 cm
- 10-535-12WB Tungsten carbide, curved, pointed, 11.5 cm, Wavecut



Reynolds scissors

- 20-609-15B Tungsten carbide, curved, blunt, tapered tips, 16 cm
- 20-609-15WB Tungsten carbide, curved, blunt, tapered tips, 16 cm, Wavecut
- 20-610-15B Curved, blunt, tapered tips, 16 cm
- 20-610-15WB Curved, blunt, tapered tips, 16 cm, Wavecut



Metzenbaum scissors

- 20-613-15B Tungsten carbide, straight, blunt, fine, 14.5 cm
- 20-613-15WB Tungsten carbide, straight, blunt, fine, 14.5 cm, Wavecut
- 20-613-18B Tungsten carbide, straight, blunt, fine, 18.0 cm
- 20-613-18WB Tungsten carbide, straight, blunt, fine, 18.0 cm, Wavecut
- 20-613-20B Tungsten carbide, straight, blunt, fine, 20.0 cm
- 20-613-20WB Tungsten carbide, straight, blunt, fine, 20.0 cm, Wavecut
- 20-615-15B Tungsten carbide, curved, blunt, fine, 14.5 cm
- 20-615-15WB Tungsten carbide, curved, blunt, fine, 14.5 cm, Wavecut
- 20-615-18B Tungsten carbide, curved, blunt, fine, 18.0 cm
- 20-615-18WB Tungsten carbide, curved, blunt, fine, 18.0 cm, Wavecut
- 20-615-20B Tungsten carbide, curved, blunt, fine, 20.0 cm
- 20-615-20WB Tungsten carbide, curved, blunt, fine, 20.0 cm, Wavecut



Joseph scissors

- 20-602-14B Tungsten carbide, straight, pointed, 14 cm
- 20-602-14WB Tungsten carbide, straight, pointed, 14 cm, Wavecut
- 20-603-14B Tungsten carbide, curved, pointed, 14 cm
- 20-603-14WB Tungsten carbide, curved, pointed, 14 cm, Wavecut



Kilner rhinoplasty scissors

- 20-622-13B Tungsten carbide, curved, blunt, flat tip, 13.5 cm
- 20-622-13WB Tungsten carbide, curved, blunt, flat tip, 13.5 cm, Wavecut
- 20-622-15B Tungsten carbide, curved, blunt, flat tip, 15.0 cm
- 20-622-15WB Tungsten carbide, curved, blunt, flat tip, 15.0 cm, Wavecut



Stevens vascular scissors

- 20-623-11B straight, tapered tips, 11.5 cm
- 20-623-11WB straight, tapered tips, 11.5 cm, Wavecut
- 20-624-11B curved, tapered tips, 11.5 cm
- 20-624-11WB curved, tapered tips, 11.5 cm, Wavecut



Black Line

Ultimate Edition


Walter scissors

20-630-10B angled, 11.5 cm



Preyer retractor

10-446-11B slim blade, studded on both sides, detachable, 11 cm



Cottle lower lateral scissors

20-632-12B strongly curved, 12.5 cm



Fomon nasal scissors

20-635-13WB angled, fine, 14.0 cm, Wavecut

Heymann nasal scissors

20-680-16WB straight, fine, 16.0 cm, Wavecut

20-683-18WB straight, robust, 18.0 cm, Wavecut

Blakesley nasal forceps

20-215-00B straight, 3.0 mm, 12 cm

20-215-01B straight, 3.5 mm, 12 cm

20-215-02B straight, 4.0 mm, 12 cm

20-215-03B straight, 5.0 mm, 12 cm

20-215-04B straight, 6.0 mm, 12 cm

20-225-00B 45° curved upwards, 3.0 mm, 12 cm

20-225-01B 45° curved upwards, 3.5 mm, 12 cm

20-225-02B 45° curved upwards, 4.0 mm, 12 cm

20-225-03B 45° curved upwards, 5.0 mm, 12 cm

20-235-01B 90° curved upwards, 3.5 mm, 12 cm

20-235-02B 90° curved upwards, 4.0 mm, 12 cm


20-235-03B 90° curved upwards, 5.0 mm, 12 cm

20-225-11B 45° curved upwards, 3.5 mm, 15 cm

20-225-12B 45° curved upwards, 4.0 mm, 15 cm

20-225-14B 45° curved upwards, 6.0 mm, 15 cm





© 2019 by SPIGGLE & THEIS Medizintechnik GmbH
Modifications and errors - including pictures - reserved.
Please note also that not all goods might be for sale
in your country.



Burghof 14 · 51491 Overath · Germany
Phone: +49 2206 90 81 - 0
Fax: +49 2206 90 81 - 13
info@spiggle-theis.com
www.spiggle-theis.com