

Inspiring for Splendid Healthcare

ISO | ISO | CE
9001:2015 | 13485:2016



OSSEOUS PRODUCT CATALOG

- OTOLOGICAL MIDDLE EAR PROSTHESIS
- OTOLOGICAL VENTILATION TUBES
- RHINOLOGY PRODUCTS
- MEDICAL PVA DRESSINGS
- LARYNGOLOGY PRODUCTS
- ENT INSTRUMENTS
- ENT PLASMA SYSTEM - Coaglator II
- ABSORBABLE GELATIN SPONGES



MediTop
— Medical Products —



Material: CP Titanium Medical Grade ASTM F 67

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E1104	1.40	2.55	2.55	1.20	10 Units, Gamma Sterile
E1105	1.30	2.55	2.55	1.20	10 Units, Gamma Sterile
E1106	1.50	2.75	2.75	1.20	10 Units, Gamma Sterile

Material: CP Titanium Medical Grade ASTM F 67

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E1151	1.14	2.40	2.40	2.20	10 Units, Gamma Sterile



Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2109	1.14	3.54	2.40	3.85	10 Units, ETO Sterile
E2109B	1.14	3.54	2.40	3.85	10 Units, ETO Sterile
E2109T	1.14	3.54	2.40	3.85	10 Units, ETO Sterile



Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2110	1.14	3.46	2.50	3.04	10 Units, ETO Sterile
E2110B	1.14	3.46	2.50	3.04	10 Units, ETO Sterile
E2110T	1.14	3.46	2.50	3.04	10 Units, ETO Sterile



Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2114	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114B	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114T	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114W	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114BW	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114BT	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2138	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138W	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138T	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138B	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138BW	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138BT	1.00	2.4	2.4	2.2	10 Units, ETO Sterile



Shepard & Mini Shepard

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2137	1.14	2.3	2.3	1.6	10 Units, ETO Sterile
E2137B	1.14	2.3	2.3	1.6	10 Units, ETO Sterile



Reuter Bobbin, Reuter Bobbin Blue

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2136	1.27	2.8	2.8	2.2	10 Units, ETO Sterile
E2136B	1.27	2.8	2.8	2.2	10 Units, ETO Sterile



Baxter Beveled Bobbin, Baxter Beveled Bobbin Blue

B - Blue, T - Tail, W - Wire, All dimensions are in MM

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2112	1.14	3.4	2.3	2.4	10 Units, ETO Sterile
E2112B	1.14	3.4	2.3	2.4	10 Units, ETO Sterile
E2112T	1.14	3.4	2.3	2.4	10 Units, ETO Sterile
E2112W	1.14	3.4	2.3	2.4	10 Units, ETO Sterile
E2112BW	1.14	3.4	2.3	2.4	10 Units, ETO Sterile
E2112BT	1.14	3.4	2.3	2.4	10 Units, ETO Sterile
E2113	0.76	2.2	1.5	1.5	10 Units, ETO Sterile
E2113B	0.76	2.2	1.5	1.5	10 Units, ETO Sterile
E2113T	0.76	2.2	1.5	1.5	10 Units, ETO Sterile
E2113W	0.76	2.2	1.5	1.5	10 Units, ETO Sterile



Shah & Mini Shah



Donaldson,
Donaldson Blue

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2124	1.14	2.3	2.3	2.2	10 Units, ETO Sterile
E2124B	1.14	2.3	2.3	2.2	10 Units, ETO Sterile

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2117	1.14	2.8	2.8	1.7	10 Units, ETO Sterile
E2117B	1.14	2.8	2.8	1.7	10 Units, ETO Sterile



Bevel Bobbin,
Bevel Bobbin Blue



Straight Tube

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2116	1.14	2.7	1.54	7.25	10 Units, ETO Sterile

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2123	1.27	3.0	3.0	2.15	10 Units, ETO Sterile
E2123B	1.27	3.0	3.0	2.15	10 Units, ETO Sterile
E2123T	1.27	3.0	3.0	2.15	10 Units, ETO Sterile
E2125	1.14	2.50	2.50	2.0	10 Units, ETO Sterile
E2125B	1.14	2.50	2.50	2.0	10 Units, ETO Sterile
E2125T	1.14	2.50	2.50	2.0	10 Units, ETO Sterile



Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E3111	1.10	2.3	2.3	2.2	10 Units, ETO Sterile

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E3126	1.10	3.5	2.5	1.65	10 Units, ETO Sterile



Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Total Length	Units/Box
E3115S	1.14	8.0	6.0	10 Units, ETO Sterile
E3115M	1.14	8.0	9.0	10 Units, ETO Sterile
E3115L	1.14	8.0	12.0	10 Units, ETO Sterile

Material : Silver Plated Gold

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E1170	1.25	2.55	2.55	1.20	10 Units, Gamma Sterile



B - Blue, T - Tail, W - Wire, All dimensions are in MM

EON MEDITECH

 **MediTop**
— Medical Products —

 **EON MEDITECH** Inspiring for Splendid Healthcare®
PVT. LTD.

www.eonmeditech.com



EON MEDITECH PVT. LTD.
12, Gurukrupa Society, Utran Road,
Surat - 394105, Gujrat, INDIA.
Ph.: +91 261-2494330
E-mail : info@eonmeditech.com



CMC MEDICAL DEVICES & DRUGS S.L.
C/Horacio Lengo No. 18, CP 29006,
Malaga, Spain.
E-mail : info@cmcmedicaldevices.com

