

Premium Scissors

The next generation



Quality you can rely on

30
Years
SPIGGLE & THEIS
1994 - 2024



Setting new standards

Premium Scissors


Scissors are an indispensable element of modern surgery. They represent one of the most important instruments in daily OR use. Countless designs exist to meet daily requirements.

Scissors **MADE IN GERMANY** represent superlative quality, precision, and durability.

With our **Premium Line**, we are setting new standards that revolutionize cutting behavior *in a way never before imagined*.

Experience a new generation of scissors whose precision, quality and durability cannot be beaten. Make your own impression and contact us.

We look forward to hearing from you.




Premium scissors

Excellent sharpness

» In many years as an ENT surgeon, I've never worked with **better** and **sharper** scissors than with the new **SPIGGLE & THEIS Premium Scissors**. «

Prof. Livio Presutti, University Hospital Bologna


Precision, sharpness and optimal **haptic qualities** are the characteristics for our new generation of Premium Scissors. Thanks to their special design and surface texture, our Premium Scissors allow excellent **cutting strength**. Bruised or torn tissue is a thing of the past.



Advantages:

- Extreme sharpness
- Excellent haptic qualities
- Smooth cutting

Quality
Made in
Germany



High gloss polished



Blade grinding



Metzenbaum scissors

blunt-blunt



Metzenbaum Scissors, Premium Line

Art.-No.	Length (cm)	Description
20-630-14	14.5	straight
20-630-18	18.0	straight
20-630-20	20.0	straight
20-631-14	14.5	curved
20-631-18	18.0	curved
20-631-20	20.0	curved
20-633-14	14.5	fine, straight
20-633-18	18.0	fine, straight
20-633-20	20.0	fine, straight
20-634-14	14.5	fine, curved
20-634-18	18.0	fine, curved
20-634-20	20.0	fine, curved

Metzenbaum-Baby scissors

blunt-blunt



Metzenbaum-Baby Scissors, Premium Line

Art.-No.	Length (cm)	Description
20-646-12	11.5	straight
20-642-12	11.5	curved

Reynolds scissors

curved, blunt-blunt, tapered tip



Reynolds Scissors, Premium Line

Art.-No.	Length (cm)
20-635-15	15.0
20-635-18	18.0

Joseph scissors
pointed-pointed



Joseph Scissors, Premium Line

Art.-No.	Length (cm)	Description
20-640-14	14.5	straight
20-641-14	14.5	curved

Scissors
pointed-pointed



Scissors, Premium Line

Art.-No.	Length (cm)	Description
20-647-12	11.5	straight
20-648-12	11.5	curved

Heymann nasal scissors
serrated, straight, fine



Heymann nasal scissors, Premium Line

Art.-No.	Length (cm)	Description
20-650-16	16.0	fine
20-651-18	18.0	thick



